

SALL2 antibody

Product Information

Catalog No.:	FNab07582
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Probable transcription factor that plays a role in eye development before, during, and after optic fissure closure.

Immunogen information

Immunogen:	sal-like 2(Drosophila)
Synonyms:	Sal-like protein 2 Zinc finger protein 795 Zinc finger protein SALL2 Zinc finger protein Spalt-2 (Sal-2, hSal2) SALL2 KIAA0360 SAL2 ZNF795
Observed MW:	105 kDa
Uniprot ID :	Q9Y467

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab07582(SALL2 antibody) at dilution of 1:300