

## SAMD4B antibody

### Product Information

Catalog No.:	FNab07589
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Has transcriptional repressor activity. Overexpression inhibits the transcriptional activities of AP-1, p53/TP53 and CDKN1A.

### Immunogen information

Immunogen:	sterile alpha motif domain containing 4B
Synonyms:	Protein Smaug homolog 2 (Smaug 2, hSmaug2) Sterile alpha motif domain-containing protein 4B (SAM domain-containing protein 4B) SAMD4B SMAUG2
Observed MW:	70-75 kDa
Uniprot ID :	Q5PRF9

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab07589(SAMD4B antibody) at dilution of 1:800