

SAMD9L antibody

Product Information

Catalog No.:	FNab07591
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	sterile alpha motif domain containing 9-like
Synonyms:	Sterile alpha motif domain-containing protein 9-like (SAM domain-containing protein 9-like) SAMD9L C7orf6 DRIF2 KIAA2005 UEF
Observed MW:	160-180 kDa
Uniprot ID :	Q8IVG5

Application

Reactivity:	Human
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab07591(SAMD9L Antibody) at dilution of 1:50



COLO 320 cells were subjected to SDS PAGE followed by western blot with FNab07591(SAMD9L Antibody) at dilution of 1:1000