

SAP130 antibody

Product Information

Catalog No.:	FNab07595
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Acts as a transcriptional repressor. May function in the assembly and/or enzymatic activity of the mSin3A corepressor complex or in mediating interactions between the complex and other regulatory complexes.

Immunogen information


Immunogen:	Sin3A-associated protein, 130kDa
Synonyms:	Histone deacetylase complex subunit SAP130 130 kDa Sin3-associated polypeptide Sin3-associated polypeptide p130 SAP130
Observed MW:	130 kDa
Uniprot ID :	Q9H0E3

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas cancer using FNab07595(SAP130 antibody) at dilution of 1:50



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab07595(SAP130 antibody) at dilution of 1:400