

MAPK13 antibody

Product Information

Catalog No.:	FNab07597
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	mitogen-activated protein kinase 13
Synonyms:	Mitogen-activated protein kinase 13 (MAP kinase 13, MAPK 13) Mitogen-activated protein kinase p38 delta (MAP kinase p38 delta) Stress-activated protein kinase 4 MAPK13 PRKM13 SAPK4
Observed MW:	42 kDa
Uniprot ID :	O15264

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab07597(SAPK4 antibody) at dilution of 1:300