

SAR1B antibody

Product Information

Catalog No.:	FNab07603
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Involved in transport from the endoplasmic reticulum to the Golgi apparatus. Activated by the guanine nucleotide exchange factor PREB. Involved in the selection of the protein cargo and the assembly of the COPII coat complex.

Immunogen information


Immunogen:	SAR1 homolog B
Synonyms:	GTP-binding protein SAR1b GTP-binding protein B (GTBPB) SAR1B SARA2 SARB
Observed MW:	25-27 kDa
Uniprot ID :	Q9Y6B6

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer using FNab07603(SAR1B antibody) at dilution of 1:100



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab07603(SAR1B antibody) at dilution of 1:500