

ZFYVE9 antibody

Product Information

Catalog No.:	FNab07605
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	zinc finger, FYVE domain containing 9
Synonyms:	Zinc finger FYVE domain-containing protein 9 Mothers against decapentaplegic homolog-interacting protein (Madh-interacting protein) Novel serine protease (NSP) Receptor activation anchor (hSARA) Smad anchor for receptor activation ZFYVE9 MADHIP SARA SMADIP
Observed MW:	156 kDa, 180 kDa
Uniprot ID :	O95405

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF, FC, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IF: 1:20-1:200
Image:	



IP Result of anti-SARA (IP:FNab07605, 5ug;
Detection:FNab07605 1:500) with mouse brain
tissue lysate 4000ug.



mouse brain tissue were subjected to SDS PAGE
followed by western blot with
FNab07605(ZFYVE9 antibody) at dilution of
1:1000