

## SBNO1 antibody

### Product Information

Catalog No.:	FNab07618
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	strawberry notch homolog 1(Drosophila)
Synonyms:	Protein strawberry notch homolog 1 Monocyte protein 3 (MOP-3) SBNO1 MOP3
Observed MW:	145-154 kDa
Uniprot ID :	A3KN83

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:200
Image:	



Immunohistochemistry of paraffin-embedded mouse testis using FNab07618(SBNO1 antibody) at dilution of 1:50



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab07618( SBNO1 Antibody) at dilution of 1:600