

P3H4 antibody

Product Information

Catalog No.:	FNab07619
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	synaptonemal complex protein SC65
Synonyms:	Endoplasmic reticulum protein SC65 Leprecan-like protein 4 Nucleolar autoantigen No55 Prolyl 3-hydroxylase family member 4 Synaptonemal complex protein SC65 P3H4 LEPREL4 NOL55 SC65
Observed MW:	55 kDa
Uniprot ID :	Q92791

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



human bladder tissue were subjected to SDS
PAGE followed by western blot with
FNab07619(SC65 antibody) at dilution of 1:400