

SCML1 antibody

Product Information

Catalog No.: FNab07639

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months(Avoid repeated freeze / thaw cycles.)

Background

Putative Polycomb group(PcG) protein. PcG proteins act by forming multiprotein complexes, which are required to maintain the transcriptionally repressive state of homeotic genes throughout development. May be involved in spermatogenesis during sexual maturation(By similarity).

Immunogen information

Immunogen: sex comb on midleg-like 1(Drosophila)

Synonyms: Sex comb on midleg-like protein 1|SCML1

Observed MW: 23 kDa Uniprot ID: Q9UN30


Application

Reactivity: Human, Mouse Tested Application: ELISA, WB, IF

Recommended dilution: WB: 1:500-1:2000; IF: 1:20-1:200

Image:





HepG2 cells were subjected to SDS PAGE followed by western blot with FNab07639(SCML1 antibody) at dilution of 1:400