

## SCRN1 antibody

### Product Information

Catalog No.:	FNab07657
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Regulates exocytosis in mast cells. Increases both the extent of secretion and the sensitivity of mast cells to stimulation with calcium(By similarity).

### Immunogen information

Immunogen:	secernin 1
Synonyms:	Secernin-1 SCRN1 KIAA0193
Observed MW:	46 kDa
Uniprot ID :	Q12765

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:200-1:1000
Image:	



A549 cells were subjected to SDS PAGE followed by western blot with FNab07657(SCRN1 Antibody) at dilution of 1:300