

CWC27 antibody

Product Information

Catalog No.:	FNab07664
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	serologically defined colon cancer antigen 10
Synonyms:	Spliceosome-associated protein CWC27 homolog Antigen NY-CO-10 Probable inactive peptidyl-prolyl cis-trans isomerase CWC27 homolog (PPIase CWC27) Serologically defined colon cancer antigen 10 CWC27 SDCCAG10
Observed MW:	54 kDa
Uniprot ID :	Q6UX04

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab07664(SDCCAG10 antibody) at dilution of 1:500