

SDSL antibody

Product Information

Catalog No.:	FNab07675
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Has low serine dehydratase and threonine dehydratase activity.

Immunogen information


Immunogen:	serine dehydratase-like
Synonyms:	Serine dehydratase-like L-serine deaminase L-serine dehydratase/L-threonine deaminase L-threonine dehydratase (TDH) Serine dehydratase 2 (SDH 2) SDSL
Observed MW:	55 kDa
Uniprot ID :	Q96GA7

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver using FNab07675(SDSL antibody) at dilution of 1:100



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab07675(SDSL antibody) at dilution of 1:700