

## SEC22B antibody

### Product Information

Catalog No.:	FNab07679
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


SNARE involved in targeting and fusion of ER-derived transport vesicles with the Golgi complex as well as Golgi-derived retrograde transport vesicles with the ER.

### Immunogen information

Immunogen:	SEC22 vesicle trafficking protein homolog B
Synonyms:	Vesicle-trafficking protein SEC22b ER-Golgi SNARE of 24 kDa (ERS-24, ERS24) SEC22 vesicle-trafficking protein homolog B SEC22 vesicle-trafficking protein-like 1 SEC22B SEC22L1
Observed MW:	24 kDa
Uniprot ID :	O75396

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab07679(SEC22B antibody) at dilution of 1:500