

## SEC61B-Specific antibody

### Product Information

Catalog No.:	FNab07691
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	Sec61 beta subunit
Synonyms:	Protein transport protein Sec61 subunit beta SEC61B
Observed MW:	10 kDa
Uniprot ID :	P60468

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:1000-1:4000; IF: 1:20-1:200; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human cervical cancer tissue slide using FNab07691(SEC61B-Specific Antibody) at dilution of 1:50



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab07691(SEC61B-Specific antibody) at dilution of 1:500