

SCG5 antibody

Product Information

Catalog No.:	FNab07700
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	secretogranin V(7B2 protein)
Synonyms:	Neuroendocrine protein 7B2 Pituitary polypeptide Secretogranin V Secretogranin-5 Secretory granule endocrine protein I N-terminal peptide C-terminal peptide SCG5 SGNE1
Observed MW:	27 kDa
Uniprot ID :	P05408

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas cancer using FNab07700(SCG5 antibody) at dilution of 1:100



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab07700(SCG5 antibody) at dilution of 1:2000