

## SECTM1 antibody

### Product Information

Catalog No.:	FNab07701
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	3H2
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May be involved in thymocyte signaling.

### Immunogen information


Immunogen:	secreted and transmembrane 1
Synonyms:	Secreted and transmembrane protein 1 Protein K-12 SECTM1 K12
Observed MW:	27 kDa
Uniprot ID :	Q8WVN6

### Application

Reactivity:	Human
Tested Application:	ELISA, WB,IHC, FC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using FNab07701(SECTM1 Antibody) at dilution of 1:200



K-562 cells were subjected to SDS PAGE followed by western blot with FNab07701(SECTM1 Antibody) at dilution of 1:1000