

DHCR24 antibody

Product Information

Catalog No.:	FNab07703
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	24-dehydrocholesterol reductase
Synonyms:	Delta(24)-sterol reductase 24-dehydrocholesterol reductase 3-beta-hydroxysterol Delta-24-reductase Diminuto/dwarf1 homolog Seladin-1 DHCR24 KIAA0018
Observed MW:	60 kDa
Uniprot ID :	Q15392

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer using FNab07703(DHCR24 antibody) at dilution of 1:100



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab07703(DHCR24 antibody) at dilution of 1:800