

---

## SELENON antibody

### Product Information

Catalog No.:	FNab07716
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	selenoprotein N, 1
Synonyms:	Selenoprotein N (SeIN) SELENON SELN SEPN1
Observed MW:	70 kDa
Uniprot ID :	Q9NZV5

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human heart tissue slide using FNab07716(SEPN1 Antibody) at dilution of 1:50



A549 cells were subjected to SDS PAGE followed by western blot with FNab07716(SEPN1 Antibody) at dilution of 1:300