

MSRB1 antibody

Product Information

Catalog No.:	FNab07730
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	selenoprotein X, 1
Synonyms:	Methionine-R-sulfoxide reductase B1 (MsrB1) Selenoprotein X (SelX) MSRB1 SEPX1
Observed MW:	None
Uniprot ID :	Q9NZV6

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemical of paraffin-embedded mouse kidney using FNab07730(SEPX1 antibody) at dilution of 1:50