
SERINC2-Specific antibody

Product Information

Catalog No.:	FNab07736
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	serine incorporator 2
Synonyms:	Serine incorporator 2 Tumor differentially expressed protein 2-like SERINC2 TDE2L
Observed MW:	65 kDa
Uniprot ID :	Q96SA4

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:1000
Image:	



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab07736(SERINC2-Specific antibody) at dilution of 1:100