

## SERINC3 antibody

### Product Information

Catalog No.:	FNab07737
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Restriction factor required to restrict infectivity of lentiviruses, such as HIV-1: acts by inhibiting an early step of viral infection. Impairs the penetration of the viral particle into the cytoplasm(PubMed:26416733, PubMed:26416734).

### Immunogen information


Immunogen:	serine incorporator 3
Synonyms:	Serine incorporator 3 Tumor differentially expressed protein 1 SERINC3 DIFF33 TDE1
Observed MW:	67 kDa
Uniprot ID :	Q13530

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lung cancer using FNab07737(SERINC3 antibody) at dilution of 1:50



mouse pancreas tissue were subjected to SDS PAGE followed by western blot with FNab07737(SERINC3 antibody) at dilution of 1:500