

Serpina3k antibody

Product Information

| | |
|---------------|--|
| Catalog No.: | FNab07741 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

Background


Contrapsin inhibits trypsin-like proteases.

Immunogen information

| | |
|--------------|--|
| Immunogen: | serine(or cysteine) peptidase inhibitor, clade A, member 3K |
| Synonyms: | Serine protease inhibitor A3K (Serpina A3K) Contrapsin SPI-2 Serpina3k Mcm2 Spi2 |
| Observed MW: | 50-60 kDa |
| Uniprot ID : | P07759 |

Application

| | |
|-----------------------|------------------|
| Reactivity: | Mouse |
| Tested Application: | ELISA, WB |
| Recommended dilution: | WB: 1:500-1:2000 |
| Image: | |



mouse heart tissue were subjected to SDS PAGE followed by western blot with FNab07741(Serpina3k antibody) at dilution of 1:1000