

SERPINB9 antibody

Product Information

Catalog No.:	FNab07749
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Granzyme B inhibitor.

Immunogen information

Immunogen:	serpin peptidase inhibitor, clade B(ovalbumin), member 9
Synonyms:	Serpin B9 Cytoplasmic antiproteinase 3 (CAP-3, CAP3) Peptidase inhibitor 9 (PI-9) SERPINB9 PI9
Observed MW:	42 kDa
Uniprot ID :	P50453

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



K-562 cells were subjected to SDS PAGE followed by western blot with FNab07749(SERPINB9 antibody) at dilution of 1:2000