

## SERPINE2 antibody

### Product Information

Catalog No.:	FNab07752
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	4C3
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Serine protease inhibitor with activity toward thrombin, trypsin, and urokinase. Promotes neurite extension by inhibiting thrombin. Binds heparin.

### Immunogen information

Immunogen:	serpin peptidase inhibitor, clade E(nexin, plasminogen activator inhibitor type 1), member 2
Synonyms:	Glia-derived nexin (GDN) Peptidase inhibitor 7 (PI-7) Protease nexin 1 (PN-1) Protease nexin I Serpine E2 SERPINE2 PI7 PN1
Observed MW:	44 kDa
Uniprot ID :	P07093

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab07752(SERPINE2 Antibody) at dilution of 1:1000