

## SF3B5 antibody

### Product Information

Catalog No.:	FNab07781
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	splicing factor 3b, subunit 5, 10kDa
Synonyms:	Splicing factor 3B subunit 5 (SF3b5) Pre-mRNA-splicing factor SF3b 10 kDa subunit SF3B5 SF3B10
Observed MW:	10-12 kDa
Uniprot ID :	Q9BWJ5

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



human adrenal gland tissue were subjected to SDS PAGE followed by western blot with FNab07781(SF3B5 antibody) at dilution of 1:400