

AKAP17A antibody

Product Information

Catalog No.:	FNab07789
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	splicing factor, arginine/serine-rich 17A
Synonyms:	A-kinase anchor protein 17A (AKAP-17A) 721P B-lymphocyte antigen Protein XE7 Protein kinase A-anchoring protein 17A (PRKA17A) Splicing factor, arginine/serine-rich 17A AKAP17A CXYorf3 DXYS155E SFRS17A XE7
Observed MW:	81 kDa
Uniprot ID :	Q02040

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas cancer using FNab07789(SFRS17A antibody) at dilution of 1:100



human kidney tissue were subjected to SDS PAGE followed by western blot with FNab07789(SFRS17A antibody) at dilution of 1:400