

SRSF7 antibody

Product Information

Catalog No.:	FNab07792
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	splicing factor, arginine/serine-rich 7, 35kDa
Synonyms:	Serine/arginine-rich splicing factor 7 Splicing factor 9G8 Splicing factor, arginine/serine-rich 7 SRSF7 SFRS7
Observed MW:	35 kDa
Uniprot ID :	Q16629

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, IF, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human prostate cancer using FNab07792(SFRS7 antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab07792(SFRS7 Antibody) at dilution of 1:1000