

SRSF9 antibody

Product Information

Catalog No.:	FNab07794
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	splicing factor, arginine/serine-rich 9
Synonyms:	Serine/arginine-rich splicing factor 9 Pre-mRNA-splicing factor SRp30C Splicing factor, arginine/serine-rich 9 SRSF9 SFRS9 SRP30C
Observed MW:	26 kDa
Uniprot ID :	Q13242

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab07794(SFRS9 antibody) at dilution of 1:1000