

## ARHGEF26 antibody

### Product Information

Catalog No.:	FNab07804
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	Src homology 3 domain-containing guanine nucleotide exchange factor
Synonyms:	Rho guanine nucleotide exchange factor 26 SH3 domain-containing guanine exchange factor ARHGEF26 SGEF
Observed MW:	97 kDa
Uniprot ID :	Q96DR7

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab07804(SGEF antibody) at dilution of 1:300