

## SGK3 antibody

### Product Information


Catalog No.:	FNab07808
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information


Immunogen:	serum/glucocorticoid regulated kinase family, member 3
Synonyms:	Serine/threonine-protein kinase Sgk3 Cytokine-independent survival kinase Serum/glucocorticoid-regulated kinase 3 Serum/glucocorticoid-regulated kinase-like SGK3 CISK SGKL
Observed MW:	57 kDa
Uniprot ID :	Q96BR1

### Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas cancer using FNab07808(SGK3 antibody) at dilution of 1:100



IP Result of anti-SGK3 (IP:FNab07808, 4ug; Detection:FNab07808 1:500) with mouse heart tissue lysate 3200ug.



HeLa cells were subjected to SDS PAGE followed by western blot with FNab07808(SGK3 antibody) at dilution of 1:300