

SLC5A2 antibody

Product Information

Catalog No.:	FNab07809
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	solute carrier family 5(sodium/glucose cotransporter), member 2
Synonyms:	Sodium/glucose cotransporter 2 (Na(+)/glucose cotransporter 2) Low affinity sodium-glucose cotransporter Solute carrier family 5 member 2 SLC5A2 SGLT2
Observed MW:	69-70 kDa
Uniprot ID :	P31639

Application

Reactivity:	Human
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney tissue slide using FNab07809(SLC5A2 Antibody) at dilution of 1:200



A549 cells were subjected to SDS PAGE followed by western blot with FNab07809(SLC5A2 Antibody) at dilution of 1:600