

SH3BGR antibody

Product Information

Catalog No.:	FNab07826
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	SH3 domain binding glutamic acid-rich protein
Synonyms:	SH3 domain-binding glutamic acid-rich protein (SH3BGR protein) 21-glutamic acid-rich protein (21-GARP) SH3BGR
Observed MW:	26-30 kDa
Uniprot ID :	P55822

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human heart tissue slide using FNab07826(SH3BGR Antibody) at dilution of 1:50



mouse heart tissue were subjected to SDS PAGE followed by western blot with FNab07826(SH3BGR Antibody) at dilution of 1:1000