

SH3GL2 antibody

Product Information

Catalog No.:	FNab07833
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Implicated in synaptic vesicle endocytosis. May recruit other proteins to membranes with high curvature.

Immunogen information


Immunogen:	SH3-domain GRB2-like 2
Synonyms:	Endophilin-A1 EEN-B1 Endophilin-1 SH3 domain protein 2A SH3 domain-containing GRB2-like protein 2 SH3GL2 CNSA2 SH3D2A
Observed MW:	40 kDa
Uniprot ID :	Q99962

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human gliomas using FNab07833(SH3GL2 antibody) at dilution of 1:50



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab07833(SH3GL2 Antibody) at dilution of 1:1000