

SHISA4 antibody

Product Information

Catalog No.:	FNab07848
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	shisa homolog 4(Xenopus laevis)
Synonyms:	Protein shisa-4 Transmembrane protein 58 SHISA4 C1orf40 TMEM58
Observed MW:	25-27 kDa
Uniprot ID :	Q96DD7

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab07848(SHISA4 antibody) at dilution of 1:2000