

SHKBP1 antibody

Product Information

Catalog No.:	FNab07849
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Inhibits CBL-SH3KBP1 complex mediated down-regulation of EGFR signaling by sequestration of SH3KBP1. Binds to SH3KBP1 and prevents its interaction with CBL and inhibits translocation of SH3KBP1 to EGFR containing vesicles upon EGF stimulation.

Immunogen information

Immunogen:	SH3KBP1 binding protein 1
Synonyms:	SH3KBP1-binding protein 1 SETA-binding protein 1 SHKBP1 SB1
Observed MW:	76 kDa
Uniprot ID :	Q8TBC3

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



A431 cells were subjected to SDS PAGE followed by western blot with FNab07849(SHKBP1 antibody) at dilution of 1:300