

SHOX antibody

Product Information

Catalog No.:	FNab07852
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Controls fundamental aspects of growth and development.

Immunogen information

Immunogen:	short stature homeobox
Synonyms:	Short stature homeobox protein Pseudoautosomal homeobox-containing osteogenic protein Short stature homeobox-containing protein SHOX PHOG
Observed MW:	25-32 kDa
Uniprot ID :	O15266

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab07852(SHOX Antibody) at dilution of 1:1000