

GEMIN2 antibody

Product Information

Catalog No.:	FNab07871
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	survival of motor neuron protein interacting protein 1
Synonyms:	Gem-associated protein 2 (Gemin-2) Component of gems 2 Survival of motor neuron protein-interacting protein 1 (SMN-interacting protein 1) GEMIN2 SIP1
Observed MW:	32-35 kDa
Uniprot ID :	O14893

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab07871(SIP1 antibody) at dilution of 1:500