

SIPA1 antibody

Product Information

Catalog No.:	FNab07872
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


GTPase activator for the nuclear Ras-related regulatory proteins Rap1 and Rap2 in vitro, converting it to the putatively inactive GDP-bound state.

Immunogen information

Immunogen:	signal-induced proliferation-associated 1
Synonyms:	Signal-induced proliferation-associated protein 1 (Sipa-1) GTPase-activating protein Spa-1 p130 SPA-1 SIPA1 SPA1
Observed MW:	130 kDa
Uniprot ID :	Q96FS4

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



Raji cells were subjected to SDS PAGE followed by western blot with FNab07872(SIPA1 antibody) at dilution of 1:3000