

SIPA1L1 antibody

Product Information

Catalog No.:	FNab07873
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Stimulates the GTPase activity of RAP2A. Promotes reorganization of the actin cytoskeleton and recruits DLG4 to F-actin. Contributes to the regulation of dendritic spine morphogenesis(By similarity).

Immunogen information

Immunogen:	signal-induced proliferation-associated 1 like 1
Synonyms:	Signal-induced proliferation-associated 1-like protein 1 (SIPA1-like protein 1) High-risk human papilloma viruses E6 oncoproteins targeted protein 1 (E6-targeted protein 1) SIPA1L1 E6TP1 KIAA0440
Observed MW:	200 kDa
Uniprot ID :	O43166

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab07873(SIPA1L1 Antibody) at dilution of 1:1500