

SIRT2-Specific antibody

Product Information

Catalog No.:	FNab07880
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	sirtuin(silent mating type information regulation 2 homolog) 2
Synonyms:	NAD-dependent protein deacetylase sirtuin-2 NAD-dependent protein defatty-acylase sirtuin-2 Regulatory protein SIR2 homolog 2 SIR2-like protein 2 SIRT2 SIR2L SIR2L2
Observed MW:	43 kDa
Uniprot ID :	Q8IXJ6

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB, IP
Recommended dilution:	WB: 1:500-1:1000; IP: 1:500-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human heart tissue slide using FNab07880(SIRT2 Antibody) at dilution of 1:50



IP result of anti-SIRT2-Specific(IP:FNab07880, 4ug; Detection:FNab07880 1:800) with mouse brain tissue lysate 4000 ug



rat brain tissue were subjected to SDS PAGE followed by western blot with FNab07880(SIRT2 antibody) at dilution of 1:600