

SLAMF1 antibody

Product Information

Catalog No.:	FNab07902
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	signaling lymphocytic activation molecule family member 1
Synonyms:	Signaling lymphocytic activation molecule CDw150 IPO-3 SLAM family member 1 SLAMF1 SLAM
Observed MW:	70 kDa
Uniprot ID :	Q13291

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:200-1:1000
Image:	



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab07902(SLAM/CD150 antibody at dilution of 1:300