

SLC22A23 antibody

Product Information

Catalog No.:	FNab07924
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	solute carrier family 22, member 23
Synonyms:	Solute carrier family 22 member 23 SLC22A23 C6orf85
Observed MW:	66-74 kDa
Uniprot ID :	A1A5C7

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



L02 cells were subjected to SDS PAGE followed by western blot with FNab07924(SLC22A23 Antibody) at dilution of 1:1500