

SLC25A10 antibody

Product Information

Catalog No.:	FNab07928
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Involved in translocation of malonate, malate and succinate in exchange for phosphate, sulfate, sulfite or thiosulfate across mitochondrial inner membrane.

Immunogen information


Immunogen:	solute carrier family 25(mitochondrial carrier; dicarboxylate transporter), member 10
Synonyms:	Mitochondrial dicarboxylate carrier (DIC) Solute carrier family 25 member 10 SLC25A10 DIC
Observed MW:	32 kDa
Uniprot ID :	Q9UBX3

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:1000; IP: 1:500-1:1000
Image:	



IP Result of anti-SLC25A10 (IP: FNab07928, 3ug;
Detection: FNab07928 1:800) with mouse liver
tissue lysate 3000ug.



mouse liver tissue were subjected to SDS PAGE
followed by western blot with
FNab07928(SLC25A10 antibody) at dilution of
1:600