

## SLC25A15 antibody

### Product Information

Catalog No.:	FNab07931
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Ornithine transport across inner mitochondrial membrane, from the cytoplasm to the matrix.

### Immunogen information


Immunogen:	solute carrier family 25(mitochondrial carrier; ornithine transporter) member 15
Synonyms:	Mitochondrial ornithine transporter 1 Solute carrier family 25 member 15 SLC25A15 ORC1 ORNT1
Observed MW:	30-33 kDa
Uniprot ID :	Q9Y619

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human heart using FNab07931(SLC25A15 antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab07931(SLC25A15 antibody) at dilution of 1:500