

SLC25A46 antibody

Product Information

Catalog No.:	FNab07944
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May play a role in mitochondrial dynamics by controlling mitochondrial membrane fission.

Immunogen information


Immunogen:	solute carrier family 25, member 46
Synonyms:	Mitochondrial outer membrane protein SLC25A46 Solute carrier family 25 member 46 SLC25A46 TB1
Observed MW:	46 kDa
Uniprot ID :	Q96AG3

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human medulloblastoma using FNab07944(SLC25A46 antibody) at dilution of 1:50



MCF-7 cells were subjected to SDS PAGE followed by western blot with FNab07944(SLC25A46 Antibody) at dilution of 1:600