

## SLC27A3 antibody

### Product Information

Catalog No.:	FNab07949
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Has acyl-CoA ligase activity for long-chain and very-long-chain fatty acids. Does not exhibit fatty acid transport activity(By similarity).

### Immunogen information


Immunogen:	solute carrier family 27(fatty acid transporter), member 3
Synonyms:	Long-chain fatty acid transport protein 3 (FATP-3, Fatty acid transport protein 3) Arachidonate--CoA ligase Long-chain-fatty-acid--CoA ligase Solute carrier family 27 member 3 Very long-chain acyl-CoA synthetase homolog 3 (VLCS-3) SLC27A3 ACSVL3 FATP3
Observed MW:	70-79 kDa
Uniprot ID :	Q5K4L6

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:2000
Image:	



IP Result of anti-SLC27A3 (IP:FNab07949, 4ug;  
Detection:FNab07949 1:500) with mouse testis  
tissue lysate 4000ug.



human testis tissue were subjected to SDS PAGE  
followed by western blot with  
FNab07949(SLC27A3 antibody) at dilution of  
1:1500