

## SLC5A3 antibody

### Product Information

Catalog No.:	FNab07971
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Prevents intracellular accumulation of high concentrations of myo-inositol(an osmolyte) that result in impairment of cellular function.

### Immunogen information

Immunogen:	solute carrier family 5(sodium/myo-inositol cotransporter), member 3
Synonyms:	Sodium/myo-inositol cotransporter (Na(+)/myo-inositol cotransporter) Sodium/myo-inositol transporter 1 (SMIT1) Solute carrier family 5 member 3 SLC5A3
Observed MW:	77 kDa
Uniprot ID :	P53794

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab07972(SLC5A3 antibody) at dilution of 1:300