

## SLFNL1 antibody

### Product Information

Catalog No.:	FNab07980
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	schlafen-like 1
Synonyms:	Schlafen-like protein 1 SLFNL1
Observed MW:	46 kDa
Uniprot ID :	Q499Z3

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney using FNab07980(SLFNL1 antibody) at dilution of 1:100



COLO 320 cells were subjected to SDS PAGE followed by western blot with FNab07980(SLFNL1 antibody) at dilution of 1:800