

SLK antibody

Product Information

Catalog No.:	FNab07985
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Mediates apoptosis and actin stress fiber dissolution.

Immunogen information

Immunogen:	STE20-like kinase(yeast)
Synonyms:	STE20-like serine/threonine-protein kinase (STE20-like kinase, hSLK) CTCL tumor antigen se20-9 STE20-related serine/threonine-protein kinase (STE20-related kinase) Serine/threonine-protein kinase 2 SLK KIAA0204 STK2
Observed MW:	95 kDa
Uniprot ID :	Q9H2G2

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab07985(SLK antibody) at dilution of 1:2000